

Body fat scale manual

One. OKOK APP Downloading

Downloading body finger scale APP There are several ways:

1. scan two-dimensional code download



Available on the
App Store



OKOK Healthcare

Two. APP install



Three. Open APP registration information

1. The registration information can be completed in the following ways

2. 1. login the number of the mobile phone number
3. 2. Log in with QQ, WeChat, mailbox account



Connect the phone

- 1、 Turn on Bluetooth on your phone

When using the OKOK software, you need to turn on the scale first, then connect via Bluetooth Click on the "OKOK app downloaded on the phone," and the following interface will appear, then You can step onto the body fat scale, and once the scale locks the current weight, the mobile app will If the following interface appears, click "Confirm Binding".



APP Function Interface

After measurement, the data will be automatically synchronized



Specification

Display : All data displayed when power on

Weighing range: 0.2~180kg

Division; 1D=100g/0.2lb

Unit; KG/LB/ST

Voltage: 4.5V

Low battery indication: 2.4V to 0.2V, display "Lo" and power off

Current: ≤25mA

Standby current: ≤10uA

Overload indication: overload more than 2.5kgs, total 182.5kgs display "Err"

Best working temperature: 5~35℃

Storage temperature: -10~60℃

Moisture Environment: 90%RH

Automatic shutdown time: 8 seconds

Automatic zero point tracking

Power supply: If the purchased battery model is 2*1.5V (AA batteries)

If the purchased charging model is USB-powered at 3.7V

WARRANTY CARD

Customer information

Product name		Product model	
Date of purchase		Shop purchase	
Dealer Seal			
Customer name		Customer phone	
Customer address			

certificate

Model: P-01

Name :human health scale

examination clerk: 08

